

**PHILIPS**


**VOLCANO**

*ChromaFlo*


Imaging

Enhanced clarity.  
**Simple.**


Easy assessment of stent apposition, lumen size, and more.  
ChromaFlo imaging highlights blood flow red at the touch of a button.




Lumen sizing




Bifurcation




Malapposed stent



Thrombus




Bioresorbable scaffold



Dissection

# ChromaFlo is available on Philips Volcano Core integrated, Core Mobile, s5, and s5i systems

## Chroma.



## Simple to use

- Simply press the button
- No special procedures, no flushing during imaging
- Works with Philips Volcano's fast, plug-and-play digital IVUS catheters as shown below

## Broad versatility

- Appropriate for coronary and peripheral vessels including left mains, bifurcations, SFA and iliac arteries
- Provides both lumen and vessel information
  - Colorization highlights the lumen
  - Automatically co-registered grayscale IVUS allows observation of plaque burden, calcium, and vessel size

## Compatible IVUS catheters



### Eagle Eye Platinum ST catheter

0.014" guide wire compatible



### Eagle Eye Platinum catheter

0.014" guide wire compatible



### Visions PV .014P catheter

0.014" guide wire compatible



### Visions PV .018 catheter

0.018" guide wire compatible

